

## Oil Red O Stain, Propylene Glycol - Technical Memo

### SOLUTION:

Oil Red O Stain, Propylene Glycol

250 ml

Part 12772A

500 ml

Part 12772B

### Additionally Needed:

Formalin 10%, Phosphate Buffered

Part 1090

Propylene Glycol 100%, ACS

Part 13391

Propylene Glycol 85%, Aqueous

Part 133912

Hematoxylin Stain, Mayer Modified

Part 1202

Lithium Carbonate, Saturated Aqueous

Part 12215

or Scott Tap Water Substitute

Part 1380

Mount-Quick Aqueous

Part 6271A

*For storage requirements and expiration date refer to individual bottle labels.*

### APPLICATION:

Newcomer Supply Oil Red O Stain, Propylene Glycol procedure is classified as a physical staining method and is used for identification of fat/lipid in frozen sections.

### METHOD:

**Fixation:** Fresh tissue or formalin fixed unprocessed tissue

a. See Procedure Note #1.

**Technique:** Frozen sections cut at 8-10 microns mounted on adhesive slides

**Solutions:** All solutions are manufactured by Newcomer Supply, Inc.

All Newcomer Supply stain procedures are designed to be used with Coplin jars filled to 40 ml following the staining procedure provided below.

### STAINING PROCEDURE:

- Fix frozen section slides in Formalin 10%, Phosphate Buffered (1090) for 1 minute.
  - See Procedure Note #2.
- Rinse sections carefully in two changes of distilled water.
- Blot off excess water and place slides in Propylene Glycol 100%, ACS (13391) for 2-5 minutes.
- Place slides directly into Oil Red O Stain, Propylene Glycol for 1 hour. Agitate occasionally or place Coplin jar on rotator/shaker with continuous gentle agitation.
  - See Procedure Notes #3 and #4.
- Differentiate in Propylene Glycol 85%, Aqueous (133912) with agitation for a minimum of 3 minutes.
- Rinse gently in two changes of distilled water.
- Counterstain with Hematoxylin Stain, Mayer Modified (1202) for 2-3 minutes, depending on preference of nuclear stain intensity.
- Wash gently in several changes of tap water.
- Blue slides in Lithium Carbonate, Saturated Aqueous (12215) or Scott Tap Water Substitute (1380) for 10 dips.
  - See Procedure Note #5.
- Wash gently in several changes of tap water.
- Blot excess water from slide; coverslip with Mount-Quick Aqueous (6271A) mounting medium.
  - See Procedure Note #6.

### RESULTS:

Fat      Bright red  
Nuclei    Blue to dark blue

### PROCEDURE NOTES:

- To freeze formalin fixed unprocessed tissue post-fixation: Place specimen in tissue cassette, wash in running tap water for 5 minutes. Remove tissue from cassette and blot well removing all excess water from tissue. Freeze tissue (fresh or formalin fixed) according to laboratory protocol.
- Frozen formalin fixed tissue does not require an additional formalin fixation step.
- To decrease staining time; preheat Oil Red O Stain, Propylene Glycol in a 60°C oven and decrease incubation time to 7-10 minutes.
- If a filmy precipitate develops in Oil Red O Stain, Propylene Glycol, filter with coarse filter paper.
- The use of a bluing agent in this procedure is an optional step.
- Use minimal pressure when applying coverslip or fat/lipid staining may be disturbed. To remove trapped air bubbles under coverslip; soak slide in warm water until coverslip can easily be removed. Blot excess water from slide and remount with new coverslip and Mount-Quick Aqueous mounting medium.

### REFERENCES:

- Carson, Freida L., and Christa Hladik. *Histotechnology: A Self-Instructional Text*. 3rd ed. Chicago, Ill.: American Society of Clinical Pathologists, 2009. 184-186
- Prophet, Edna B., Bob Mills, Jacquelyn Arrington, and Leslie Sobin. *Laboratory Methods in Histotechnology*. Washington, D.C.: American Registry of Pathology. 1992.178.
- Sheehan, Dezna C., and Barbara B. Hrapchak. *Theory and Practice of Histotechnology*. 2nd ed. St. Louis: Mosby, 1980. 205.
- Modifications developed by Newcomer Supply Laboratory.