



ANIMAL ON ANIMAL IHC



Tired of needing multiple staining kits, an array of normal sera and blocking agents for each animal tissue?

STOP THE MADNESS!!!

(and experience the joy of background free animal IHC)

JUST THREE INCREDIBLE PRODUCTS FOR ALL YOUR ANIMAL TISSUES!
ALL SPECIES AND ALL ORGANS

1 - Stat-Q Staining Kit:

- The kit of choice for all animal IHC!
- Multispecies 3-step staining kit.
- Produces strong, clean results.



2 - Background Buster:

- Get rid of background staining easily.
- Makes child's play of Mouse on Mouse IHC!
- No more cocktails of normal sera & blocking agents.

3 - Enhancing Wash Buffer:

(Optional to perk up your staining intensity)

- Use in place of your current buffer. (no additional steps)
- Gives a turbo boost to your IHCs!
- Eliminate or minimize high heat epitope retrieval.



Timings:

Background Buster

10 minutes (*Human Tissue*)
or
30 minutes (*Mouse on Mouse*)

Primary Ab

30 minutes

Link

10 minutes

Label

10 minutes

Chromogen

5 minutes

ALSO EXCELLENT FOR HUMAN TISSUE, ALL DAY & EVERY DAY!

STAT-Q STAINING KIT AND BACKGROUND BUSTER ARE A MUST FOR BACKGROUND FREE ANIMAL IHC!

Description	Volume	Slide Capacity	Price
Stat-Q™ Peroxidase Multispecies Staining System (DAB or AEC) <i>(Other Monovalent links available for Stat-Q Staining Kit)</i>	20 ml	200-300	\$345.00
	60 ml	800-1,000	\$645.00
Stat-Q™ Alk Phos Multispecies Staining System (Fast Red or Innovex Brown) <i>(Other Monovalent links available for Stat-Q Staining Kit)</i>	20 ml	200-300	\$345.00
	60 ml	800-1,000	\$645.00
Background Buster	50ml	700-1,000	\$175.00
	125ml	1,400-2,000	\$275.00
Enhancing Wash Buffer - HRP (Amplifying AEC & DAB stains)	1 liter		\$89.00
	4 liter		\$195.00
Enhancing Wash Buffer - Alk Phos. (Amplifying Fast Red & Innovex Brown Stains)	1 liter		\$89.00
	4 liter		\$195.00

Distributed by:



2505 Parview Road
Middleton, WI 53562

Phone (800) 383-7799
Fax (608) 831-0866
www.newcomersupply.com

A Product of

