

BACKGROUND STAINING ON YOUR IMMUNO SLIDES JUST MET ITS MATCH...

MEET THE BACKGROUND BUSTER BOYZ!

When your slides become a mess, just give us a call and we'll do the rest!

When you have finally had enough of non-specific background staining on your IHC slides, give us a call and we'll introduce you to the *Buster Boyz*.

These two are feared by any background staining!

- ✓ From primary antibodies
- ✓ From Staining Kits
- ✓ From Substrate/Chromogens

(Just the mere mentioning of the Buster Boyz in the presence of your immuno slides might start to clear them up.)



THE POWER OF BACKGROUND BUSTER



A powerful recombinant protein technology for eliminating general background staining!



Allows staining of identical species antibodies and tissue (i.e. mouse Ab on mouse tissue).



Replaces the use of all normal serum, powdered milk, casein and other blocking agents.



Excellent for both frozen and paraffin sections.

A MUST for Animal Tissue Staining!!!

<u>Description</u>	<u>Volume</u>	<u>Slide Capacity</u>	<u>Order #</u>	<u>Price</u>
Background Buster (Shelf life of 3 years)	50ml	700-1,000 slides	NB306-50	\$175.00
	125ml	1,400-2,000 slides	NB306	\$275.00



800-383-7799

Distributed by:

NEWCOMER SUPPLY™

2505 Parview Road
Middleton, WI 53562

Phone (800) 383-7799
Fax (608) 831-0866
www.newcomersupply.com

A Product of

INNOVEX
BIOSCIENCES