

**T
H
E

C
H
R
O
M
O
G
E
N

E
N
H
A
N
C
E
R
S
!**

**GET THE MOST OUT OF YOUR CHROMOGENS,
EVERYTIME!**

You will be the envy of all your co-workers when your immunos always look their best!



Take control of those sometimes "difficult" antibodies and enhance their staining intensity with an easy incubation step!

It's as simple as:

- 1) Incubate with the chromogen of choice
- 2) Rinse
- 3) Apply the Chromogen Enhancer and incubate for 5 to 10 minutes
- 4) Rinse
- 5) Counterstain and mount!

The Enhancers will:

- Intensify & Amplify your stain
- Minimize false negative staining
- Insure your chromogens are at their best!

The Enhancers will NOT:

- Create any background or non-specific staining

Think of it as an insurance policy for your peace of mind.



800-383-7799

<u>Description</u>	<u>Volume</u>	<u>Shelf Life</u>	<u>Order #</u>	<u>Price</u>
AEC Chromogen Enhancer	15ml	3 years	NB319	\$135.00
Fast Red Enhancer	15ml	3 years	NB322	\$155.00
Quick DAB Enhancer	60ml	3 years	NB308	\$195.00

Distributed by:



2505 Parview Road
Middleton, WI 53562

Phone (800) 383-7799
Fax (608) 831-0866
www.newcomersupply.com

A Product of

