

TAKE CONTROL OF YOUR LYMPHOID TISSUE

Did you know those pesky little white blood cells in your tissue sections play with your immuno results? Yep, it's those darn Fc receptors that bind to the Fc region of the primary and secondary antibodies. The background and false positives they produce are especially troublesome when you're staining CD markers, lymphomas, bone marrow and blood smears.

Your Frustration is Over:

INTRODUCING

FC RECEPTOR BLOCKER!

Also Works on
Brain Tissue!



- ☺ One 20 minute incubation at room temp. It's that Easy!
- ☺ Eliminates false positive staining in lymphomas, CD marker stainings, bone marrow, lymph nodes & cytosmears.

Put a Smile on your pathologist's face with
Easy to Read Kappa and Lambda staining!



Fc Receptor Block is a MUST for accurate lymphoma and leukemia typing.

<u>Description</u>	<u>Volume</u>	<u>Slide Capacity</u>	<u>Order #</u>	<u>Price</u>
Fc Receptor Blocker (Shelf life of 3 years)	15ml	150-250 slides	NB309-15	\$189.00
	60ml	700-1,000 slides	NB309	\$485.00



800-383-7799

Distributed by:

NEWCOMER SUPPLY™

2505 Parview Road
Middleton, WI 53562

Phone (800) 383-7799
Fax (608) 831-0866
www.newcomersupply.com

A Product of

**INNOVEX
BIOSCIENCES**