




THE CONVENIENCE OF QUICKLOAD™ SLEEVES FOR YOUR TISSUE CASSETTE PRINTER



LOAD THEM BY HAND, WHY???

...You purchased that labeler for faster through-put and convenience, it only makes sense to automate the entire process!!!



Slimsette™
Processing/Embedding
Part #5121

Slimsette™ Processing / Embedding Cassettes

- * Hinged lid that opens from the back of the cassette for easy access
- * 45° angle writing surface
- * Recessed cover to save stacking space
- * 114 openings each measuring 1 x 5mm

Slimmer size brings an abundance of advantages: (total height 6mm vs. 7mm for most other cassettes)

- Keeps grossed tissue 5mm maximum in thickness - Better Fixation!
- Allows for more cassettes to be run in the processor - more per stack!
- Allows for more cassettes to be stacked in the automatic labeler - more efficient!



Slimsette™
Biopsy
Part# 5122

Slimsette™ Biopsy / Embedding Cassettes

- * Hinged lid that opens from the back of the cassette for easy access
- * 45° angle writing surface
- * Recessed cover to save stacking space
- * 420 openings each measuring 1 x 1mm



Micromesh™
Biopsy
Part# 5123

Micromesh™ Biopsy / Embedding Cassettes

- * Hinged lid that opens from the back of the cassette for easy access
- * 45° angle writing surface
- * Recessed cover to save space
- * 1676 openings each measuring 0.38mm
- * Square compartment measures 27 x 27mm

<u>Description</u>	<u>Volume</u>	<u>Order #</u>	<u>Price</u>
Slimsette™ Processing/Embedding Cassettes in QuickLoad™ Sleeves (Multiple colors available, 10 sleeves of 75 cassettes/case)	750 / cs	5121	\$250.00!
Slimsette™ Biopsy/Embedding Cassettes in QuickLoad™ Sleeves (Multiple colors available, 10 sleeves of 75 cassettes/case)	750 / cs	5122	\$250.00!
Micromesh™ Biopsy/Embedding Cassettes in QuickLoad™ Sleeves (Multiple colors available, 10 sleeves of 75 cassettes/case)	750 / cs	5123	\$250.00!

Wow!

Distributed by:



2505 Parview Road
Middleton, WI 53562

Phone (800) 383-7799

Fax (608) 831-0866

www.newcomersupply.com