

ADHESIVE MICROSCOPE SLIDE SAMPLE FORM



Name:
Company:
Address:

City: State: ZIP:
Phone: ext:
Email:

- Samples ship to U.S. and Canada only
- Return the sample form by:
email: dprine@newcomersupply.com
fax: 608-831-0866

NewSilane

Part #: 5070

Adhesion: Very Strong

Application: Difficult and routine IHC.
Difficult tissues with long incubations.

Unit: 72 slides/box, 1,440 slides/case

Price: \$59.70/box, \$898.00/case

MAS

Part #: 5077

Adhesion: Very Strong

Surface: Hydrophilic

Application: Excellent for the most
challenging tissues and applications.

Unit: 100 slides/box, 1,000 slides/case

Price: \$295.00/case

TOMO

Part #: 5081

Adhesion: Strong

Surface: Hydrophilic

Application: For automated IHC platforms such
as Ventana/Roche Benchmark and Leica Bond.

Unit: 100 slides/box, 1,000 slides/case

Price: \$280.00/case

MAS-GP

Part #: 5079

Adhesion: Strong

Surface: Hydrophilic

Application: For all IHC applications, HIER,
pH6 / pH9 protocols, routine histology, and
frozen sections.

Unit: 100 slides/box, 1,000 slides/case

Price: \$280.00/case

APS

Part #: 5078

Adhesion: Strong

Surface: Hydrophobic, treated to provide a
permanent positive charge.

Application: IHC and routine histology.

Unit: 100 slides/box, 1,000 slides/case

Price: \$180.00/case

TRUBOND 380

Part #: 5080

Adhesion: Very Strong

Surface: Hydrophilic

Application: For HIER processed when tissue
type requires extra strong adhesion, epitope
enhancement and DNA probe procedures.

Unit: 100 slides/box, 1,000 slides/case

Price: \$280.00/case