

**FREE
SAMPLES!**

Adhesive Microscope Slides

Check boxes for the samples
you would like

APS

- Surface is treated to provide a permanent positive charge surface
- Recommended for: All IHC applications, routine histology
- Surface: Hydrophobic
- Adhesion: Strong
- \$180.00/case (1,000 slides/case)



TRUBOND 380

- Useful in preventing tissue loss under high temp. and high pH conditions
- Recommended for: HIER processes when tissue type requires extra strong adhesion, epitope enhancement & DNA probe procedures, autopsy & brain sectioning
- Surface: Hydrophilic
- Adhesion: Very Strong
- \$280.00/case (1,000 slides/case)



NEWSILANE

- Useful for "tricky" routine and IHC applications
- Difficult tissues combined with long incubations
- Extended antigen retrieval
- Available in White
- 1-10 boxes \$57.00/box; 11-19 boxes \$48.90/box; 20+ boxes \$42.90/box (approx. 72/box)



MAS-GP (Formerly TRUBOND 360)

- Useful in preventing tissue loss under high temperature and high pH conditions
- Recommended for: All IHC applications, HIER pH6 and pH9 protocols, routine histology, frozen sectioning
- Surface: Hydrophilic
- Adhesion: Strong
- \$280.00/case (1,000 slides/case)



TOMO

- Useful in preventing tissue loss under high temperature and high pH conditions
- Recommended for: Automated IHC platforms (Ventana Benchmark, Leica Bond, etc.)
- Surface: Hydrophilic
- Adhesion: Strong
- \$280.00/case (1,000 slides/case)



CREST

- Strong adhesion without background
- Surface: Hydrophilic
- \$280.00/case (10 boxes/case)

MAS

- Excellent for the most challenging tissues and applications
- Surface: Hydrophilic
- \$40.00/box (100 slides); \$280.00/case (10 boxes/case)

DIAMOND WHITE

- Produced from superior quality white glass providing clarity, low self fluorescence and uniform thickness.
- Ideal for use in histology, cytology and pathology departments
- Dimensions: 25x75mm x 1 mm thick
- Available in assorted colors
- \$365.00/case (1,440 slides/case)



SHIP SAMPLES TO:

Name _____

Company _____

Address _____

Address _____

City, ST Zip _____, _____

Phone (_____) _____ - _____ ext. _____

Email Address _____

FAX this completed form to: **608-831-0866**

OR

Email to: **dprine@newcomersupply.com**

Website