


# Which Eosin is best for you?



	<u>1082 Eosin Phloxine Set</u>	<u>1070 1% Alcoholic</u>	<u>1080 1% Aqueous</u>	<u>1072 Working Eosin</u>
<b>Dilution</b>	100ml 1% Eosin Stock 10ml 1% Phloxine B 780ml 95% ETOH	100ml stock to 300ml 95% ETOH	100ml stock to 300ml 95% ETOH	No Dilution Necessary
<b>Acetic Acid</b>	Add 4ml	.5ml/100ml	.5ml/100ml	Not required
<b>Recommended Use</b>	Phloxine adds more red to cytoplasmic staining	Cytoplasmic stain for a less intense counter- stain than 1% aqueous	Cytoplasmic stain for a more intense counter- stain than 1% alcoholic	Easy to use no dilution based on 1% aqueous formula
<b>Shelf Life</b>	1 year	2 years	1 year	1 year

